3D cell culture, assay automation and 3D printing demonstrated with a novel xeno-free hydrogel derived from wood



Introduction

GrowDex[®] is wood-based nanofibrillar cellulose (NFC) hydrogel developed for 3D cell culture. It is biocompatible with human cells and tissues but as a plant based product it does not contain any animal or human derived material.

GrowDex efficiently supports 3D cell _ physically resembling growth by extracellular matrix (ECM). The structure and mechanical properties of -GrowDex can be tuned to fulfill the requirements of different cell types and it allows the diffusion nutrients and oxygen. The hydrogel can be completely degraded to soluble glucose by enzyme treatment while retaining the 3D cellular structure. Under stress GrowDex has shear thinning properties, which make it a pipettable ready-to-use hydrogel.

Key properties of GrowDex include:

- **Ready-to-use**
- **Xeno-free material**
- **Biocompatible**
- Adjustable stiffness
- **Enzyme degradable Non-autofluoresecent No batch variation**

Stem Cells

The natural stem cell niche is a cells and the shifting of 3D culture to dynamic 3D environment supporting 2D platforms for various downstream stem cell proliferation. GrowDex applications, such as supports the proliferation of human embryonic stem cells (hESC) and human induced pluripotent stem cells karyotypes of the stem cells, shown (**hiPSC**) without feeder cells. Stem cells by chromosomal G-band analyses. form 3D spheroids when cultured in 0.5% GrowDex hydrogel (Fig. 5) and the pluripotency is maintained during the 3D culture. Enzymatic degradation mesenchymal stem cells (hMSC) to of the hydrogel with cellulase enzyme enables simple sub-culturing of the

directed differentiation. 3D cell culturing and enzymatic degradation do not affect

GrowDex has also been shown to support 3D growth human of from spheroids and establish cell function (Fig. 6).



Figure 1. GrowDex hydrogel.



Figure 5. A) GrowDex supports proliferation and 3D growth of human induced pluripotent stem cells. B) Live/Dead analysis shows the viability of the cells in spheroids. Images from Yan-Ru Lou (Ph.D.), University of Helsinki, Finland.



Figure 6. Bone marrow-derived mesenchymal stem cells form 3D spheroids in GrowDex in 24h.

Hepatic Cells

Human hepatic cell lines, such as Human primary hepatocytes are used HepaRG and HepG2 proliferate and e.g. in *in vitro* models for studying the form multicellular 3D spheroids in liver toxicity of new drug compounds. GrowDex hydrogel (Figs.2 and 3A-B). Human primary hepatocytes form 3D Cells remain viable in GrowDex, spheroids in GrowDex during 28 days indicated by Live/Dead staining of the cultures (Fig.2B). The differentiation in HepG2 spheroids can be observed as *in* toxicity studies. vivo -like cell polarization based on the concentrated presence of F-actin (Fig.

culture (Fig.3C). Long culture period enables long term drug exposure and



Removal of GrowDex by Cellulase Enzyme

100 um

Growdex can be completely removed As the cellular structure remains and complex at 37°C until the hydrogel has completely degraded.

by cellulase enzyme. UPM cellulase intact detailed downstream analysis enzyme is a purified and optimized of the formed spheroids or cellular enzyme that specifically degrades structure is possible. In addition to cellulose to soluble glucose without brightfield and confocal imaging, affecting the cells. The use of enzyme immunohistochemistry and scanning enables easy degradation of the electron microscopy (SEM) are also matrix whilst retaining the 3D possible following enzyme treatment. **spheroid structure**. The enzyme can be SEM images show that after enzyme simply mixed with cell culture medium treatment microvilli structures typical incubated with cell-hydrogel for hepatocytes are revealed (Fig. 6).

3A). Various stains can be used to visualize the spheroids in GrowDex (Fig.3B HepaRG).



Figure 2. A) Phase contrast microscope image of HepaRG spheroids cultured in GrowDex. B) The viability of HepaRG cells analysed by Live/Dead staining. Images from Liisa Kanninen (Ph.D.) and Melina Malinen (Ph.D.), University of Helsinki, Finland.



Figure 3. A) Staining of nuclei and filamentous actin (F-actin) in HepG2 spheroids. B) FDA/Draq5[™] stained HepaRG spheroids. C) Human primary hepatocyte spheroid in GrowDex after 28 days culture, staining by DAPI/Phalloidin. Images A and B from Liisa Kanninen (Ph.D.) and Melina Malinen (Ph.D.), University of Helsinki, Finland.

Automated Dispensing and 3D Printing

properties GrowDex can be used in a variety of **automatic dispensing** systems applications (Fig.4B). The hydrogel has high viscosity and yield stress at rest, but starts to flow when shear force is applied, e.g. pushing hydrogel through a pipette tip. Initial high viscosity is reestablished immediately after the shear

As a consequence of its shear-thinning Fluctuations in temperature do not affect the performance of GrowDex. Thus it can be dispensed at room (Fig.4A) that are for HTS and 3D printing temperature or 37°C without issue. This characteristic makes GrowDex an ideal support matrix for cell-based 3D high throughput screening assays.

> The ability to **3D print** GrowDex opens up new possibilities in the Biomedical field, e.g.in tissue engineering.



Figure 6. Scanning electron microscopy (SEM) images of HepG2 spheroids after cellulase enzyme treatment and silica bioreplication reveals the microvilli structure typical for hepatocytes. Images from Liisa Kanninen (Ph.D.) and Yan-Ru Lou (Ph.D.), University of Helsinki, Finland.

Conclusion and Comments

The inherent properties of GrowDex make it an ideal tool for researchers working in a number of fields, from development of 3D cell models through to 3D printing of NFC hydrogel structures as demonstrated here. In summary:

GrowDex is **biocompatible** with cells and tissues

- Proliferation and 3D spheroid formation of e.g. stem cells, hepatic cells and cancer cells
- Xeno-free, the composition is clearly defined with no batch variation

force has been removed.





Figure 4. A) GrowDex can be dispensed with automated systems, such as the Biomek NXP, demonstrating its suitability for use in 3D cell-based HTS applications. B) 3D printed tubular structures created with GrowDex. The diameter of the tube is 5mm.

GrowDex is **simple to use**

- Ready to use matrix tunable with e.g. water of culture medium
- Can be used and stored at room temperature
- No crosslinking step needed
- Shear-thinning property enables automated pipetting

GrowDex can be removed with enzyme whilst retaining 3D structure Efficient recovery of spheroids and cellular structures with well-preserved shape



www.growdex.com

growdex@upm.com

Bhattacharya M., Malinen MM et al., 2012, Nanofibrillar cellulose hydrogel promotes three-dimensional liver cell culture, Journal of Controlled Release 164:3, pp. 291–298. **References:** Lou Y. et al., 2014. The use of Nanofibrillar cellulose hydrogel as a flexible three-dimensional model to culture human pluripotent stem cells, Stem cell and development. Vol. 23:4, pp. 380–392. Malinen MM, et al. 2014. Differentiation of liver progenitor cell line to functional organotypic cultures in 3D Nanofibrillar cellulose and hyaluronan-gelatin hydrogels. Biomaterials 35:5110-5121. Lou Y., Kanninen L. et al., 2015. Silica bioreplication preserves three-dimensional spheroid structures of human pluripotent stem cells and HepG2 cells, Nature Scientific Reports | 5:13635 | DOI: 10.1038/srep13635.

UPM-Kymmene Corporation, Tukholmankatu 8, Biomedicum 2U, 00290 Helsinki, Finland